


# NEOVENTURES BIOTECHNOLOGY INC.

5 YORK STREET, UNIT # 3, LONDON, ON, N6A 1A1. CANADA .519.601.7511

## Newsletter #3: Signal detection with aptamers



In this newsletter my intent is to review the various options that exist for signal detection in aptamer based diagnostics. I am going to avoid sophisticated applications such as single molecule cantilevers or channeled light waves. For now, let's stick to the basics.

It is straightforward to have a fluorophore or a biotin attached to either end of an aptamer. As a detection aptamer in sandwich system where the capture ligand is either an antibody or an aptamer both of these moieties enable direct measurement of the amount of detection aptamer bound to the target. With biotin, we apply a streptavidin/horseradish peroxidase conjugate for colour detection. The streptavidin portion of the conjugated protein binds to the biotin on the aptamer, while the HRP portion serves as the basis for the generation of colour.

Both systems work well. The use of colour is often favoured because instruments that are capable of reading colour absorption are less expensive than fluorescence readers. Sensitivity may vary depending on the matrix analyzed. It is important to remember that the colour system is an amplification of signal, as each HRP molecule bound can effect multiple colour changes.

One interesting application of fluorescently labeled aptamers is the use of fluorescence polarization. We have developed a commercial diagnostic kit for the prediction of breadmaking potential based on this approach (check out [www.risingnumber.com](http://www.risingnumber.com) for more information). We developed an aptamer that binds to a common domain across all groups of wheat storage proteins. We labeled this aptamer with a fluorophore. We combine it with the various extracted wheat protein fractions and measure the tumbling speed of the aptamer in each protein containing solution through the use of a hand held fluorescence polarization instrument. The proteins are much larger than the aptamer and thus the complex between the protein and the aptamer represents a slower tumbling speed than the aptamer alone. The beauty of this approach is that it does not require immobilization of the aptamer, reactions are performed in an entirely free system.

A key feature of aptamers is that they are single stranded oligonucleotides. This means that antisense constructs can be designed that will anneal at varying strengths depending on the degree of homology and the length of the antisense oligonucleotide. We can engineer antisense constructs with almost any binding strength for the aptamer that we want. We have found that we can obtain optimum sensitivity within detection systems when we use antisense constructs that have approximately the same binding affinity for the aptamer as the aptamer has for its target.

